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
LS-400B / LS-1250B

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The screenshot shows the 'Initial page' of a device setup interface. At the top right is a download icon. Below the title are three input fields: 'Date' with '00', 'Month' with '00', and 'Year' with '0000'. A note below the fields states: '* Please enter the date of initial use. Product warranty and A/S cannot be obtained if the input is incorrect.'

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This screenshot shows the 'Initial page' with the 'Year' field now containing '2026'. Below the input fields is a numeric keypad with digits 1 through 0, an 'Enter' button, and a backspace icon.

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This screenshot shows the 'Initial page' with a back arrow button added to the top left. The 'Year' field is '2026'. A note below the fields states: '* Please check your entry.'

The final screen displays a list of completion messages: 'Thank you for using this product.', 'Initialize the equipment.', and 'Please wait a moment.'

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Precautions Before Use

1) User

This product should be used by a person who understands the user manual or under the supervision of a responsible person.

2) AC Input Power

Before installing and using this product, please make sure to check the following items.

Confirm that the AC input power you are using matches the specifications indicated on the rear side of the equipment.

The AC input power outlet must be equipped with a grounding terminal, as shown in the figure on the right.

Please use the AC input power cable supplied with this equipment.

If it is unavoidable to use another cable, please use a cable rated at 250V, 3A or higher.

3) Maintenance and Inspection

Before performing maintenance or inspection, unplug the AC power cable from the outlet.

If inspection of the product is required, please contact the manufacturer or the place of purchase for assistance.

4) Warranty

Products purchased within one year may receive free repair service, provided that the defect is not clearly caused by consumer negligence.

Products damaged due to consumer negligence or products beyond the warranty period may be subject to paid repair service.

To receive warranty service or repair, please request service from the place of purchase or the manufacturer.

Shipping costs shall be paid separately by both the service requester and the service provider for each shipment.

5) Safety Precautions

Do not repair or modify this product without the cooperation or approval of the manufacturer.

If the product is repaired or modified without approval, the safety of the user and the performance of the product cannot be guaranteed, and the user shall be responsible for any resulting consequences.

6) Safety Symbols

The user must operate the equipment with sufficient care to avoid electric shock and must be protected from electrical hazards.

Careless use may cause damage to the equipment.

<Note> The warning and caution symbols indicated on this equipment and in the user manual are as follows.

	Warning Symbol: Electrical Hazard		
	Earth (Ground) Terminal: Functional Ground Terminal		
	Protective Conductor Terminal		
	On (Stomacher)		Off (Stomacher)

Precautions for Use

- Before use, adjust the panel nut to align it horizontally with the table.
- Do not operate the equipment under overload conditions during initial operation. Frequent warning messages may shorten the service life of the motor.
- Before bringing the paddle close to the sample, be sure to operate the LED cap. For example, start from level 1 and gradually operate the paddle so that it approaches the sample.
- Do not turn off the power while the equipment is operating.
- If an error message appears and the equipment stops, press the return button located above the LCD screen.

Table of Contents

1. Product Information

1-1. Overview	4
1-2. Basic Functions	4
1-3. Product Inspection Before Use	4
1-4. Operating Conditions	5
1-5. Cleaning Method	5
1-6. Input Power Requirements	5

2. Product Components and Description

2-1. Screen Layout and Description	6
------------------------------------	---

3. Operation

3-1. Menu Selection	8
3-2. Operation	9

4. Alarms

4-1. Door Open Alarm During Operation	11
4-2. Equipment Error Alarm	11
4-3. Motor Overload Alarm	11
4-4. LED Lighting	11

5. Functions

5-1. "Save Menu" Setting and Saving	12
5-2. "Customized Operation" Setting and Saving	15
5-3. Distance Calibration	16
5-4. Starting New Equipment	17

Product Warranty Certificate	18
------------------------------	----

1. Product Information

1-1. Overview

- The LED Stomacher LS-400B / LS-1250B finely homogenizes samples required in laboratories, allowing experiments to be performed more conveniently and accurately.
- This equipment is designed as a bench-top unit and can be used easily.
- The dimensions of the equipment are as follows.

LS - 400B	LS - 1250B
268 X 410 X 303mm, 16kg	376 X 531 X 370mm, 26kg

1-2. Basic Functions

- Timer setting function: 1 second to 2 hours
- Grinding and mixing speed setting: 1 to 10 steps, 1 to 12 strokes
- Convenient 5-inch digital LCD touch panel
- Paddle adjustment function: 30 levels, in mm units, with digital automatic paddle distance adjustment
- Program memory: Stores 9 "Save" settings and 1 "Customized" user setting according to the hardness of food samples
- 24V high-performance, low-noise motor
- Window LED: 4-color lamp status display and buzzer alarm function
- Paddles can be completely detached for cleaning

1-3. Product Inspection Before Use

- After receiving this product, inspect it for any damage that may have occurred during transportation.
- If any damage is found, please immediately contact BNF Korea Co., Ltd. (www.bnfkorea.co.kr) or the supplier of this product. We will provide replacement or repair as quickly as possible.
- A product warranty certificate is attached to the last page of this user manual. Filling in the purchase date, place of purchase, and product serial number on the warranty certificate may be useful for warranty service.
- This product undergoes several tests before shipment. During this process, minor scratches may occur on the exterior case; however, these are not subject to after-sales service.

1-4. Operating Conditions

- This equipment can be used without loss of performance when installed under the following ambient conditions.
- Ambient temperature: 0 to 40°C
- Humidity: 80% or less
- Altitude: 2,000 m or less
- Ensure sufficient clearance at the bottom and rear of the product for proper air circulation.
- Use this product in a stable environment free from vibration.

1-5. Cleaning Method


- Be sure to disconnect the power before cleaning.
- Regular cleaning is recommended to ensure long-term use of this product.
- For hygienic use, remove the paddles from the product and clean them thoroughly.
- Wipe the exterior surface with a soft cloth moistened with alcohol, detergent, or clean water.
- Do not use sharp metal objects for cleaning, as they may damage the product.
- During cleaning, protect the LCD touch panel from water and foreign substances. Failure to do so may cause malfunction or product failure.

1-6. Input Power Requirements

- This product is designed for single-phase AC 220V, 47–63Hz.
- If you need to use a power source different from the input power indicated on the rear side, be sure to contact BNF Korea Co., Ltd. (www.bnfkorea.co.kr).
- When replacing the fuse, be sure to remove the power cord from the rear panel before pulling out the fuse holder.
- Use a replacement fuse rated at 250V, 5A.

Fuse Replacement Guide

Example: For 250V 6.3A, use a 5A fuse.

Danger: Be sure to disconnect the power before changing the input power or replacing the fuse.


2. Product Components and Description

2-1. Screen Layout and Description


Initial Equipment Screen

- ① Paddle Adjustment Button
(Stabilization operation before starting the equipment)
- ② Equipment Start Button


Menu Selection Screen

- ① Saved Menu Load Button
- ② General Settings Screen Button

Saved Menu Operation Screen

- ① Speed Level Display and Adjustment
- ② Paddle Distance Level Display and Adjustment
- ③ Operation Time Display and Adjustment
- ④ Operation Time Display
- ⑤ Start Button
- ⑥ Stop Button
- ⑦ Save Settings Button
- ⑧ Saved Name Display and Setting Button
- ⑨ Button to Move to the Menu Selection Screen


Customized Operation Screen

- ① Speed Level Display and Adjustment
- ② Paddle Distance Level Display and Adjustment
- ③ Operation Time Display and Adjustment
- ④ Operation Time Display
- ⑤ Start Button
- ⑥ Stop Button
- ⑦ Save Settings Button
- ⑧ Button to Move to the Menu Selection Screen


3. Operation

3-1. Menu Selection

- The screen switches to the menu selection screen.
- The LED light of the equipment is displayed in green. If the front cover is open, it is displayed in blue.


- On the menu selection screen, you can select from **9 saved menus** and **1 customized setting**.
- You can move between screens and select a menu using the left/right navigation buttons.

User Saved Menu (Folder Icon)

- When a user saved menu is elected,
- the screen moves to the dedicated operation screen.
- The paddle distance automatically moves to the previously set and saved position.

General Operation (Customized)

- When General Operation, or Customized, is selected, the screen moves to the dedicated operation screen.
- The paddle distance automatically moves to the previously set and saved position.

3-2. Operation

*** Please check the following when installing a pack in the Stomacher.**

- ① Use a pack with an appropriate capacity.
- ② Fix the pack so that it is positioned in the center.
- ③ Turn the cover lever on the right side downward to secure the pack. The paddle strokes will start automatically.
- ④ It is recommended to first set the distance to a small value, then gradually increase it while checking the movement of the paddles.
- ⑤ Set the distance and speed appropriately to prevent overload on the equipment.

*** If motor overload is detected, an alarm will occur. In this case, adjust the distance backward and restart the operation. When restarting, adjust the distance and speed while checking the movement of the paddles.**


No.	Operation	Description
1	Speed Level Setting	  Press the left button to decrease the "Speed Level" value.  Press the right button to increase the "Speed Level" value.  The level can also be adjusted by dragging the slider. * The level can be adjusted from 1 to 10.
2	Distance Setting	  Press the left button to decrease the "Distance" value.  Press the right button to increase the "Distance" value.  The distance can also be adjusted by dragging the slider. * The distance can be adjusted from 0 to 30 mm.

3	Operation Time Setting	  Press the left button to decrease the "Time" value.  Press the right button to increase the "Time" value.  The time can also be adjusted by dragging the slider. * The time can be adjusted from 1 second to 2 hours.  The time can be adjusted by selecting hours, minutes, and seconds. Each value can be selected by clicking it, and the selected value can then be adjusted. (Hours: 1-hour increments / Minutes: 1-minute increments / Seconds: 1-second increments) The selected value is displayed in blue font.
4	Start	Fix the pack so that it is positioned in the center. Turn the cover lever on the right side downward to secure the pack. When the door is closed, the paddle strokes are automatically activated. With the door closed, press the Start button  to begin operation. The strokes will start operating. During operation, the LED light changes to white. * If an alarm occurs, the strokes stop and the LED light changes to red.
5	Operation	The equipment operates according to the set operating conditions. The distance can be adjusted during operation. The operation time can be adjusted during operation. The Speed Level can be adjusted during operation.
6	Operation Complete	Operation ends after the set time has elapsed.
7	Stop Operation	To forcibly stop the operation during a test, press the Stop button. 
8	Operation Complete Status	When operation is complete, the white LED light blinks and the buzzer sounds three times. Afterward, the light changes to green.

4. Alarms


4-1. Door Open Alarm During Operation

- If the door is opened during operation, the operation stops and an alarm occurs.
- In this case, operation can be resumed by closing the door again.
- When resuming operation, approximately 3 seconds of stabilization time is required.


4-2. Equipment Error Alarm

- If an error occurs in the equipment, an alarm will occur.
- In this case, stop using the equipment and be sure to contact the manufacturer or seller.

4-3. Motor Overload Alarm

- If an overload occurs in the motor of the equipment, an alarm will occur.
- In case ①, press the Back button (⬅️) to move to the main setting screen.
- In case ②, reset the distance backward and adjust the speed, then restart the test. Reset the distance and speed while checking the movement of the paddles.

※ If the alarm continues to occur, be sure to contact the manufacturer or seller.

4-4. LED Lighting


- The LED color of this product changes according to its status. The status indications are as follows.

Color	Status
Green	Standby before operation / Door closed
Blue	Standby before operation / Door open
White	Operating status / Door closed
Red	Door open during operation / Alarm occurs during standby before operation


5. Functions

5-1. "Save Menu" Setting and Saving

- The user can load a preset "Save Menu" and use it for operation.
- The "Save Menu" can be edited and saved on the "Save Menu Operation Screen."


No.	Operation	Description
1	Move to Menu Selection Screen	On the operation screen, press the Back button  to move to the "Menu Selection Screen." On the initial screen, press the START button  to move to the "Menu Selection Screen."
2	Select Save Menu	Select the folder icon you want to save from among the 9 options.  When the selected folder icon is pressed, the screen moves to the operation screen. The conditions previously set for that folder icon will be displayed.
3	Set Operating Conditions	Set the Speed / Distance / Time. * For the input method, refer to "3-2. Operation."

4	Set Save Name	 SAVE 01 Press the save name displayed at the upper right of the operation screen. The screen moves to the input window where the save name can be edited.   Press the uppercase/lowercase switch button to move to the uppercase/lowercase input window.   Press the number switch button to move to the number input window. 
---	------------------	--

4	Set Save Name	  Press the delete button to delete the existing folder name, then enter the desired folder name. When the folder name setting is complete, press the save button  to complete the setting. When you move to the operation screen, the folder name will be changed to the set folder name. * The folder name can be saved with up to 10 characters.
5	Save	Once all settings are complete, press the Back button  to move to the operation screen, then be sure to press the Save button  to complete the Save Menu setting. If the Save button is not clicked, the user Save Menu folder name will not be saved and the previous setting will be applied.
6	Reset	Repeat the above steps for the Save Menu you want to modify.

5-2. "Customized Operation"
Setting and Saving

- The user can save and use preset "Customized Operation" settings.
- The settings can be edited and saved on the "Customized Operation Screen."

No.	Operation	Description
1	Move to Menu Selection Screen	On the operation screen, press the Back button  to move to the "Menu Selection Screen." On the initial screen, press the START button  to move to the "Menu Selection Screen."
2	Select Customized	Select the "Customized" button.  The screen moves to the Customized Operation Screen, and the previously set conditions are displayed.
3	Set Operating Conditions	Set the Speed / Distance / Time. * For the input method, refer to "3-2. Operation."
4	Save	Once all settings are complete, press the Save button  to save the Customized Operation settings. If the Save button is not pressed, the settings will not be saved and the previous settings will be maintained.

5-3. Distance Calibration

- During initialization, the set value and the paddle position automatically move to "0 mm" and are then reset.
- Due to continuous distance adjustment, a distance error may occur in this equipment. In this case, perform the following steps for distance calibration.

No.	Operation	Description
1	Check Paddle Position and Distance Setting Value	Check the difference between the distance setting value and the actual position. If an error is detected, calibrate the distance "0" point.
2	Distance Calibration	 Press and hold the left button  of Distance for 4 seconds. When the button is pressed for more than 4 seconds, the buzzer sounds with a "beep." When the button is released, the distance moves to "0" and the value is reset. * During distance calibration, refrain from using the equipment or touching the screen, and wait until the process is complete.

5-4. Starting New Equipment

- When starting new equipment, enter the equipment purchase date.

Initial page

00

00

0000

Date

Month

Year

* Please enter the date of initial use.

Product warranty and A/S cannot be obtained if the input is incorrect.

- Enter the year, month, and day when equipment use begins, then press the Save button  to save the information.
- When you press the input field, the number input window will pop up.

Initial page

00

00

2026

Date

Month

Year

* Please enter the date of initial use.

Product warranty and A/S cannot be obtained if the input is incorrect.

1234567890

Enter

- Check the entered information.
- To make corrections, press the Back button  and edit the information.
- Once the entered information is confirmed, press the Save button  to complete the input.

←

Initial page

↓

00

00

2026

Date

Month

Year

* Please check your entry.

• Thank you for using this product.

• Initialize the equipment.

• Please wait a moment.

Product Warranty Certificate

- Warranty Period: We provide a one-year free warranty service for the mainframe from the date of purchase.
- Warranty Exclusion Policy: The following cases, including but not limited to those listed below, are not covered by the free warranty and will be subject to paid repair service.
 - a. When the warranty period has expired
 - b. Damage caused by force majeure, such as natural disasters
 - c. Damage caused by user negligence, improper use, or improper storage
 - d. Damage caused by failure to use, maintain, or adjust the product in accordance with the instructions in the user manual
 - e. Failure or damage caused by unauthorized disassembly, repair, or modification of the product; misuse; abuse; liquid ingress; accidents; or similar causes / Failure or damage caused by the use of non-genuine parts / Damage caused by damage to, forgery of, or unauthorized modification of the warranty card
 - f. Failure or damage caused by reasons other than the design, manufacturing, or quality of the product itself

Warranty Card

Product Name Product Code			
Failure Description			
Date of Purchase			
Contact Person		Contact Number	
Address			

* Note: This form is for warranty service registration. Please cut out this form and send it to the seller or manufacturer. For the detailed address of the manufacturer, please refer to the back page.

- Tel: +82-2-2666-4445
- Fax: +82-2-2666-0440
- Email: bnf@bnfkorea.net
- Website: bnfkorea.co.kr
- Online Store: bnfkorea.net
- Address: 26, Daeposandan 1-ro, Yangchon-eup, Gimpo-si, Gyeonggi-do, Republic of Korea
BNF Korea Co., Ltd.